




NEUROLOGICAL TREATMENT TABLE (NT900)



Power assisted head section



Hand held control for power assisted head section


45 years
trusted
provider .

Shown in Dusty Jade upholstery

FEATURES

The Neurological Treatment Tables have been designed for postural exercises, rolling, turning and balance training. These extra wide two section tables provide effective usage for patients who have a loss of sensation or poor balance. Also useful in treating neurological conditions; stroke, head injury and orthopaedic treatments.

In order to obtain ultimate stability and resistance to movement, we have built the NT600 & NT900 tables using lifting columns to elevate and lower the table. These synchronized lifting columns provide vertical raising and lowering which allows for total stability on the table.



NEUROLOGICAL TREATMENT TABLE

STANDARD MODEL:



2 Section Table	
Electric Hi-Lo Control	Foot Pedal and Hand Held Control
Height Adjustment	18.5" - 34" (4 lifting columns for vertical raising and lowering)
Face Cut-Out on Head Section	Keyhole cut-out with filler
Head Section Adjustment	Power assisted with hand held control
Leg Section	Fixed
Dual Hand Held Control Functionality	✓
Caster System	Dual lever retractable (5" casters)
Frame Construction	High strength seamless tubular frame
Frame Color	Durable powder coat in Gunmetal Grey or White
Upholstery Foam	1.5" thick
Antimicrobial Naugahyde Upholstery	Over 80 upholstery color options

SPECIFICATIONS:

Lifting Capacity	1000 lbs.
Dimensions	48" x 77.5" (Customizable width from 40"- 48")
Power Supply	120 Volt AC 60 Hz. Max. 5.8 Amp.
Fire Retardant Upholstery	✓

OPTIONAL ACCESSORIES:

Extra foot pedal
Extra hand held control



NEUROLOGICAL TREATMENT TABLE